Needle Sizes Animals IM SubQ IV Calf 18G x 1" 18G x 3/4" 18G x 1" 16G x 1-1/2" 14G x 2" 16G x 1" Cow 18G x 1-1/2" 16G x 1-1/2" 20G x 1-1/2" 18G x 1" Colt 20G x 1" 20G x 3/4" 18G x 1-1/2" Horse 20G x 1-1/2" 20G x 1" 20G x 1-1/2" Pig 20G x 1" 20G x 1/2" 18G x 2" Hog 18G x 1" 18G x 3/4" 18G x 4" Sow 18G x 1-1/2" 18G x 1" Sheep 18G x 1" 18G x 3/4" 18G x 1" Dog 25G x 5/8" 29G x 1/2" 25G x 5/8" Cat 25G x 5/8" 25G x 5/8" 29G x 5/8"

Improve Vaccine Results...

- Store vaccines in a cool place away from sunlight.
- Always inject vaccine into a clean, dry area.
- Reconstitute only amount of vaccine to be used within 2-3 hours.
- To avoid injusy, use only sharp needles of correct size and length.
- Inject subcutaneously (SubQ) when label permits.
- No more than 10 ml should be injected per site.
- Do not mix different products in the same syringe.
- Do not clean syringes and needles with chemical disinfectants which can inactivate vaccine.